



CHAMBRES IMPLANTABLES : TABLE DES PRODUITS



Espèce	Rats < 350gr	Rats > 350gr	Lapins	Chiens	Primates < 1-2 Kg	Primates 2 - 9 Kg	Primates > 10 Kg	Porcs	Chats	Petits animaux	Gros animaux
Taille de cathéter suggérée	- 1 french - 2 french - 3 french	- 3 french - 3,5 french - 4 french	- 5french	- 7 french - 9 french	- 3french - 3,5 french - 4 french	- 3,5 french - 4 french - 5 french	- 5 french - 7 french	- 7 french - 9 french	- 5french	- 3,5 french - 4french - 5 french	- 7 french - 9 french
Produit											
TI V-A-P				■			■	■			■
TI200V-A-P			■	■		■	■	■	■		■
Mini TI V-A-P	■	■			■	■			■	■	
GPV V-A-P			■	■		■	■	■	■		■
SLA V-A-P		■	■		■	■	■		■	■	
ROP V-A-P	■	■								■	
OmegaPort				■			■	■			■
TI-DL V-A-P				■				■			■
LCI V-A-P				■			■	■			■

Pour tout renseignement complémentaire :
 GENESTIL – 1 rue du mesnil – 60420 ROYAUCOURT
 Tel : 03 44 51 73 88
 Fax : 03 44 51 13 61
 Email : info@genestil.com

